joimax®

English

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	PAGE
1. <u>FIELD OF APPLICATION</u>	3
1.1. <u>RECOMMENDED STERILIZATION PROCEDURE</u>	3
2. <u>HANDLING</u>	4
2.1. <u>GENERAL</u>	4
2.2. <u>PREPARATIONS FOR CLEANING</u>	4
2.3. <u>FIRST USE OF A NEW STERILE CONTAINER SYSTEMS</u>	4
2.4. <u>CLEANING AND DISINFECTION</u>	4
2.4.1. <u>MANUAL CLEANING</u>	5
2.4.2. <u>MECHANICAL CLEANING</u>	5
2.4.3. <u>RECOMMENDED CLEANING PROCEDURE</u>	6
2.5. <u>INSPECTION</u>	6
2.6. <u>CHANGING THE FILTER</u>	6
3. <u>FILLING STERILE CONTAINER SYSTEMS</u>	7
4. <u>POSITIONING IN THE STERILIZER</u>	8
5. <u>FLOW CONTROL</u>	8
6. <u>STORAGE</u>	8
7. <u>MAINTENANCE/REPAIR</u>	9
8. <u>MATERIALS</u>	9
9. <u>STANDARDS APPLIED</u>	9
10. <u>WARRANTY</u>	10
11. <u>SIGNS AND SYMBOLS</u>	11

When you purchase this sterile container, you receive a high-quality product. The pages that follow explain how to use and handle it correctly. In order to keep risks and unnecessary burdens on patients, users, and third parties as low as possible, please read through the instruction manual carefully and retain it for future reference.

1. FIELD OF APPLICATION

joimax[®] sterile container systems combine proved filter technologies, tested materials, and designer characteristics in one reliable container system. They are reusable sterile container systems that come in a variety of sizes and designs in order to provide effective packaging, storage, and transportation for sterilizing instruments. The sterile container systems are perfect for fractionated vacuum procedures.

joimax[®] sterile container systems are designed to hold medical devices for sterilization. They enable and guarantee the sterilization and storage of the enclosed devices until they are used again.

The sterile container systems are available with a perforated and unperforated base, and perforated lid. The 1/1, 3/4 and 1/2 sterile container systems are also available with a safety lid.

All sterile container systems require suitable single-use paper filters which must be changed before each steam sterilization.

1.1. RECOMMENDED STERILIZATION PROCEDURE

joimax[®] sterile container systems have been tested and validated for steam sterilization with damp heat using the following sterilization parameters:

Method: 3 x prevacuum steam sterilizations
Temperature: 134 °C (273 °F)
Exposure time: 3 minutes
Drying time: 10 minutes
Contents: Standard medical instruments (scissors, clamps, pincers) and textiles

The validation covered the following joimax[®] sterile container system sizes:

- Standard model size 1/1
- Standard model size 3/4
- Standard model size 1/2

In many years of clinical experience, there are no known negative effects if the sterilization procedure is extended to the parameters below:

Method: 3 x prevacuum steam sterilizations
Temperature: 134 °C (273 °F)
Exposure time: > 18 minutes
Drying time: 20 minutes
Contents: Standard medical instruments (scissors, clamps, pincers) and textiles

2. HANDLING

2.1. GENERAL

The joimax® sterile container system comprises an aluminum alloy with an anodized surface to protect against corrosion. Aggressive cleaning agents, metal brushes, or swabs can permanently damage this surface, and therefore may not be used. If this instruction is not obeyed, the warranty will become null and void.

 **The joimax® sterile container systems can only be used by instructed or trained members of staff, in order to prevent damage to the containers, closures, seals and sterile filters.**

The joimax® sterile container systems are also available with colored lids. The colors make it easier to assign the containers to individual disciplines and specialist departments.

Sterilization indicators and colored identification labels provide information on the contents and place of use. The closure can be fitted with an anti-tamper seal which is destroyed by force on opening. Only an intact anti-tamper seal guarantees that the joimax® sterile container system was not opened by unauthorized persons.

2.2. PREPARATIONS FOR CLEANING

1. Separate the base part and lid
2. Remove contents of the sterile container system (perforated basket, instruments, etc.)
3. Remove the filter brackets from the inner side of the lid and, if necessary, the base part (for sterile container system with base perforation)
4. Dispose of the single-use filter
5. Remove the single-use anti-tamper seal and the indicator label

 **All paper filters are single-use filters and must be replaced after each insertion in the container.**

2.3. FIRST USE OF A NEW STERILE CONTAINER SYSTEM

- Before first use, the sterile container system must be thoroughly cleaned. The sterile container system must be prepared with a validated, automatic cleaning and disinfection program.
- A neutral cleaning agent (pH 7) should be used in the machine for this.
- Once preparation using the automatic cleaning and disinfection procedure is successfully completed, we recommend using steam sterilization of the devices in the fractionated steam sterilization procedure at 134 °C.
- We also recommend maintaining all moving parts on the sterile container system with instrument maintenance oil.
- After cleaning, a suitable new filter must be used (see 2.6 Changing the filter).

2.4. CLEANING AND DISINFECTION

 **Incorrect cleaning and disinfection can lead to corrosion and stress-related breakages. Therefore the manufacturer's instructions for cleaning and disinfection agents must be followed. The cleaning agents**

must be free from sodium, alkali, carbonates and must be a neutral pH 7, and mild. Only fully desalinated water (quality according to EN 285 Appendix B) is recommended for preparation of the sterile container systems.

- Before use, the sterile container system must be cleaned and disinfected.
- Containers used for disposal must be cleaned and disinfected after each use.

2.4.1. MANUAL CLEANING

- For aluminum base parts and lids, mild, neutral cleaning agents or chemical products that have been expressly approved for handling aluminum products by their manufacturing companies should be used if possible. If necessary, the devices must be checked for suitability in the appropriate procedure. For manual cleaning, a suitable damp cloth should be used. Scrubbing sponges must not be used as they destroy the surfaces and therefore passivation, and result in the warranty becoming null and void.
- After cleaning, careful rinsing with appropriate low-salt water (e.g. VE water) and sufficient drying are required.
- No metal brushes or scrubbing material may be used.

2.4.2. AUTOMATIC CLEANING

Soiling that cannot be removed regardless of the procedure in the normal cleaning process (adhesive labels, indicator strips, labels), can be removed with Eloxal cleaner. After this special procedure, the sterile container systems must be cleaned as usual. After this special procedure, the sterile container systems must be cleaned as follows:

- Neutral or other suitable cleaning and disinfection agents that have been expressly approved for preparing aluminum products. If necessary, the devices must be checked for suitability in the appropriate procedure.
- For aluminum sterile container systems, only neutral cleaning agents should be used that are expressly approved by the manufacturing companies for cleaning aluminum sterile container systems. The precise dose can also be taken from the instructions of the manufacturing company. If the program is optimized these products are also suited for cleaning surgical instruments.
- When neutralization agents are used, the products' suitability for aluminum must be checked.
- For the final rinse, desalinated water must be used.
- The cleaning devices and inserts must be suited for preparing containers and lids. This applies particularly to correctly loading inserts to ensure that the cleaning fluid can sufficiently drain unhindered and that all containers and lids dry properly.
- Reservoirs and lids may not be cleaned and disinfected while closed.
- When loading the machines, it must be ensured that fluid can flow sufficiently during the process.
- The reservoir must be positioned in the washing machine with the opening downward in order to prevent water from collecting.
- The lid must be cleaned with the inside facing down and the closures folded inward.
- After completing the automatic cleaning and disinfection process, the sterile container systems and accessories are removed with no visible residue from the fluid. However should residue be noticeable, the position of the sterile container systems and accessories in the device should be checked and changed if necessary.

2.4.3. RECOMMENDED CLEANING PROCEDURE

We recommend the following cleaning procedure, which has been validated by joimax® in a suitable cleaning and disinfection device (Miele PG 8528) with the following settings:

- 1 minute of pre-washing with cold water
- 3 minutes of cleaning with Mediclean 0.5 % (Dr. Weigert) at 45 °C
- Neutralization with demineralized water

2.5. INSPECTION

The sterile container systems must be inspected before each use for operability. Damage to the closures, seals, filter brackets, and filters, and bent and buckled parts mean that the sterile container systems must be repaired and cannot be used. Do not use damaged sterile container systems.

2.6. CHANGING THE FILTER

After the filter is changed, the filter bracket must be pressed back into the correct position until you hear a click. joimax® lids must be used with joimax® filter brackets.

- Single-use paper sterile filters must be changed before each sterilization.
- Suitability and accuracy of fit is only guaranteed by using joimax® filters.
- Warranty services can only be undertaken if original joimax® filters have been used exclusively.
- PTFE filters have been tested for a service life of 1200 cycles and must then be changed.

ATTENTION!

Only combine original joimax® individual parts such as lids, reservoirs, filters, seals and filter brackets with each other, so that impermeability and germ barriers are not jeopardized. Lids, reservoirs, filters, seals and filter brackets from other manufacturers must be approved by joimax® before they are used, otherwise joimax® makes no guarantees.

3. FILLING STERILE CONTAINER SYSTEMS

The total weight of the contents should not exceed the following recommended quantities, as otherwise satisfactory sterilization cannot be guaranteed.

Model	Dimensions L x W x H	Max. recommended contents in kg	
		Instruments	Textile load
1/1 (full) size container	580x280x100	3.8	3
	580x280x135	5.2	4.1
	580x280x150	5.8	4.6
	580x280x200	7.7	6.1
	580x280x260	10	8
3/4 size container	465x280x100	3.1	2.5
	465x280x135	4.2	3.3
	465x280x150	4.6	3.7
1/2 size container	285x280x100	1.9	1.5
	285x280x135	2.6	2
	285x280x150	2.8	2.2
	285x280x200	3.8	3
	285x280x260	4.9	3.9

When loading textiles, please make sure that the items of laundry or folded textiles are vertical. When the sterile container system is full, it should still be possible to fit your outstretched hand between the items of laundry without difficulty.

ATTENTION!

The sterilization of different sterile container systems loads and the configurations must be defined by the staff responsible for this. Endoscopes, instruments with lumen, compressed-air or mains-operated units, and instruments with canulas must be prepared for sterilization according to the manufacturer's instructions. Small baskets, boxes, or other accessories, especially with lids or flaps, should only be used together with sterile container systems if the sterile container systems are especially designed and tested for this purpose. This is confirmed by the validation of the procedure and must be complied with accordingly.

Please take further instructions on times and temperatures of sterilization from the document Instructions for Use and Preparation for Medical Devices of joimax® GmbH.

⚠ ATTENTION!

Using water-repellent insert bases (e.g. plastic/silicone mats) can lead to residual condensation inside the sterile container systems. Instead, use other suitable mats or brackets.

Check the integrity of the inserted filter and that the filter bracket closes properly. To do this, always connect the lid and reservoir using the locking mechanism. Make sure that the closure is sufficiently sealed before feeding the sterile container system into sterilization.

4. POSITIONING IN THE STERILIZER

The sterile container systems are constructed in such a way that they can be used in any commercially available large-scale sterilizer for sterilization with damp heat. Make sure that heavy sterile container systems are positioned at the bottom in the sterilization chamber. Their construction means that the sterile container systems can be easily and securely stacked on each other during sterilization without shifting around. Stacking is only recommended for sterilization cycles which run with fractionated vacuum procedures. The stacking height should not exceed 40 to 46 cm in order to achieve effective removal of air and steam penetration. The instructions of the sterilizer manufacturer must be complied with.

⚠ CAUTION! During sterilization, please note the following:

Never pack the sterile container systems into more external packaging. Never cover the perforation holes in the lid and base with any film packaging or similar, because this prevents the airflow and steam from entering the sterile container systems. As a consequence, the sterile container is distorted by a vacuum that occurs due to insufficient pressure equalization, and the sterility of the sterile container's contents cannot be guaranteed. When the sterilizer is filled and unloaded, and during transportation, the sterile container must always be carried by the handles and never the lid.

5. FLOW CONTROL

- Run the loaded sterilizer according to the instructions of the sterilizer manufacturer for the selected sterilization cycle (in relation to temperature and sterilization time). The validation results must be taken into account.
In order to minimize condensation in the sterile container, the sterile container should remain on the container carrier until it is cool enough to be handled.
- After each sterilization, the sterilization result must be assessed according to the internal instructions and validation results. The fixed connection between the lid and reservoir (closure) and the integrity of the seal made before the sterilization process must be checked before approval.

6. STORAGE

The recommended storage period for medical devices in sterilization containers should be taken from EN 58953-9 (Use of sterilization containers). Usually, the storage time depends on the storage conditions and must be defined by the responsible hygiene specialist. In the case of a particularly high requirement for asepsis, it is recommended that shorter storage periods or additional packaging are used.

Recommended storage conditions:

- Temperature: 15-26 °C
- Humidity: 30-50 %
- Air pressure: normal atmospheric pressure

The recommended storage conditions can be omitted with the series of standards.

Different container loads, storage times, and storage conditions are subject to definition by the responsible hygiene specialist. The joimax® sterile container system has been tested for a storage period of 3 months by exposure to *Bacillus subtilis* – spore suspension in an unventilated room. On this basis, we recommend a storage period of 6 weeks on open shelves and 3 months if the container is stored under protected conditions (e.g. in closed cupboards).

7. MAINTENANCE / REPAIR

- The seals have a service life of up to 500 sterilization cycles. After this the seals must be checked.
- If damage to the seals is detected, they must be immediately replaced.
- The seals should not be treated with spray, oil, or solvents. For cleaning and maintenance, it is enough to occasionally wipe them with a damp cloth.
- If damage to the sterile container systems is detected, they must be repaired.
- Maintenance and repairs to the sterile container systems may only be carried out by qualified staff. Do not try to repair the seals or fastenings yourself, so as not to damage the secure use of the containers.
- For maintenance or repairs, the sterile container systems can be sent to joimax® or an authorized repair service.
- Spare parts for 1/1, 3/4 and 1/2 sterile container systems can be purchased from joimax®:
 - Filter bracket TITF185000
 - Single-use paper filter (100 units) TITF153320

8. MATERIALS

The sterile container systems are manufactured from an anodized aluminum alloy, and the accessory from rust-free instrument steel.

9. STANDARDS APPLIED

In order to guarantee the safety of the sterile container systems during production and handling, the following standards have been taken into account:

- | | |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EN 868-2 | Packaging for terminally sterilized medical devices – Part 2: Sterilization wrap – Requirements and test methods |
| EN 868-8 | Packaging for terminally sterilized medical devices – Part 8: Re-usable sterilization containers for steam sterilizers confirming to EN 285 – Requirements and test methods |
| EN ISO 11607-1 | Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems |

EN 58952-2	Sterilization – Transport baskets for sterile barrier systems – Sterilizing baskets made of metal
EN 58952-3	Sterilization – Transport baskets for sterile barrier systems – Instrument trays for sterilizing goods made of metal
EN 58953-9	Sterilization – Sterile supply – Part 9: Use of sterilization containers
EN ISO 14937	Sterilization of health care products – General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices
ISO 11134	Sterilization of health care products – Requirements for validation and routine control – Industrial moist heat sterilization
EN ISO 17665	Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices

To guarantee sterility safety, tests were carried out by an independent and accredited testing laboratory. The goal of these tests was to validate a sterilization procedure with damp heat for the reusable joimax® sterile container systems. On the basis of the results, we recommend the sterilization procedure indicated on page 3 of these instructions for use.

10. WARRANTY

The sterile container systems were manufactured using high-quality materials and subjected to a quality test prior to delivery. Should any faults occur, contact us at the above address. The warranty shall become null and void in the event of repairs carried out by companies which are not authorized to carry out repairs by joimax®.

joimax® guarantees to manufacture the products with utmost care.

joimax® issues the legally specified warranty period with regard to the function of products. The validity period of this warranty is limited to claims made within the specified warranty period after the instrument's purchase date, with reference to repairs if required and indicating the invoice number.

This warranty only covers defects not attributable to normal wear and tear, misuse, wrong handling, no or incorrect preparation, or force majeure. This warranty excludes wear parts.

This is the only valid warranty and replaces all other issued warranty declarations.

Components of joimax® products are fully compatible. Before using individual joimax® components/products with other products, the user must ensure the user-specific compatibility of individual products.

We recommend introductory training and regular attendance at training events.

We are available for any further required information.

joimax® has no influence on the application of the product, on the diagnosis of the patient and handling of the product outside of joimax®.

joimax® cannot guarantee good results nor the complication free application of the products.

joimax® employees are not entitled to change the previously mentioned conditions, extend liability or enter additional product related obligations.

joimax® reserves the right to product changes.

This document contains information protected by copyright and property law and may not be copied either in full or in part, or transferred to a further medium in any form, without express written consent from joimax® GmbH. Forwarding to third parties is prohibited. joimax®, TES®, TESSYS®, iLESSYS®, CESSYS®, intENTS®, EndoLIF®, Percusys®, Vitegra®, Camsource®, Shrill®, Versicon®, Endovapor®, Vaporflex®, Legato®, Kyverment®, Tigrip®, Intracs® and SPOT® are registered trademarks of the company joimax®. Other products and names used here may be the registered brands of other companies. Patents pending. joimax® is grateful for any information regarding possible errors and obscurities in these instructions for use.

11. SIGNS AND SYMBOLS

	Order number
	Batch number
	Read operating instructions before use
	This CE conformity marking confirms product compliance with the Medical Devices Directive 93/42/EEC or the German Medical Devices Act (MPG).
	Manufacturer

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Copyright © 2014 joimax® GmbH. All rights reserved.
joimax®, Inc.