joimax[®]
Endoscopic Spine Experts

ENGLISH

**INSTRUCTIONS FOR USE –
RIGID ENDOSCOPES AND CERVICAL
HYBRID ENDOSCOPES**

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CHAPTER	PAGE
1. LEGAL NOTES/WARRANTY	3
2. INTRODUCTION	3
3. SAFETY NOTES/WARNING NOTES	3
4. SCOPE OF INSTRUCTIONS	6
4.1 ENDOSCOPES WITH WORKING CHANNEL	6
4.2 ENDOSCOPES WITHOUT WORKING CHANNEL	8
5. INTENDED USE	9
5.1 INTENDED USE/INTENDED PURPOSE	9
5.2 INDICATIONS AND CONTRAINDICATIONS	9
5.3 INTENDED PATIENT POPULATION	11
5.4 INTENDED USERS/USER POPULATION	11
5.5 AREAS OF APPLICATION	11
6. CONDITION OF DELIVERY	11
7. PRODUCT DESCRIPTION (ASSEMBLY/SKETCH/DESCRIPTION OF PARTS)	12
8. OPERATING OF THE DEVICES	15
9. CLEANING AND DISINFECTION	17
9.1 REPROCESSING INSTRUCTIONS	17
9.2 RESTRICTIONS AND LIMITATIONS IN REPROCESSING, DURABILITY	17
9.3 HANDLING AFTER USE	17
9.4 PREPARATION BEFORE CLEANING	18
9.5 CLEANING	19
9.5.1 AUTOMATED CLEANING/DISINFECTION (DISINFECTOR/CDD (CLEANING AND DISINFECTION DEVICE))	19
9.5.2 MANUAL CLEANING AND DISINFECTION	20
9.6 DRYING – GENERAL REQUIREMENT	21
9.7 CONTROL, MAINTANCE AND INSPECTION	21
10. STERILIZATION	23
10.1 PACKAGING	23
10.2 STEAM STERILIZATION/AUTOCCLAVING	24
10.3 PLASMA STERILIZATION	25
10.4 STORAGE	25
10.5 MATERIAL DURABILTY	25
10.6 REUSABILITY	26
10.7 ADDITIONAL INFORMATION	26
11. TROUBLESHOOTING	26
12. DISPOSAL	27
13. SERVICE AND REPAIRS	27
13.1 EMERGENCY NUMBER	28
14. SIGNS AND SYMBOLS	28

1. LEGAL NOTES / WARRANTY

Limited Warranty

joimax® guarantees to manufacture the products with utmost care.

joimax® issues the legally specified warranty period with regard to the function of instruments. The validity period of this warranty is limited to claims made within the specified warranty period after the instrument's purchase date, with reference to repairs if required and with specification of the invoice number.

This warranty only encompasses defects not caused by normal wear and tear, misuse, wrong handling, no or incorrect preparation or force majeure. This warranty excludes wearing parts. This is the only valid warranty and replaces all other issued warranty declarations.

joimax® component sets as well as individual products are compatible with each other. Before using individual joimax® products/sets with other products, the user must ensure the user-specific compatibility of individual products. We recommend introductory training and regular attendance at training events. We are available for any further required information.

joimax® has no influence on the application of the product, on the diagnosis of the patient and handling of the product outside of joimax®. joimax® cannot guarantee good results nor the complication free application of the products.

joimax® employees are not entitled to change the previously mentioned conditions, extend liability or enter additional product related obligations.

joimax® reserves the right to product changes.

This document contains information protected by copyright and property law and may not be copied in full or in parts thereof or transferred to a further medium in any form, without express written permission from joimax® GmbH. Distribution to third parties is prohibited. The products and names used here may be the registered brands of joimax® or other companies. Patents are granted or pendant.

joimax® GmbH is grateful for any information regarding possible errors and obscurities in these instructions for use.

CE-Mark

The CE marking confirms product compliance with Medical Device Regulation 2017/745.

2. INTRODUCTION

Read these instructions carefully before using your new device for the first time. This will protect you, the patient and any third parties from damages possibly caused by incorrect connection or improper operation!

Keep these instructions for use readily available.

3. SAFETY NOTES / WARNING NOTES

Incorrect use and nonobservance of these instructions, warnings and safety measures can lead to risks and severe consequences during the operation such as injury, damage or even death of the patient. The user, any third party or the endoscope itself may be harmed or injured.

These instructions for use serve to make you fully acquainted with the devices and their intended use. Therefore, they must always be readily available at the device's location.

The instructions for use contain important information required for the safe, professional and economic operation of these devices.

As future technical changes remain reserved, deviations of the contents or illustrations may be possible.

These instructions for use are intended to explain the application of endoscopes, but they do not contain instructions on the endoscopic procedure nor a detailed description of endoscopy, nor do they serve to introduce a beginner to this surgical technique!

Any person operating these devices must first read these instructions carefully.

These devices may only be operated in accordance with these instructions for use.

Use the endoscope for endoscopic purposes only!

Before each application, examine the endoscope, the respective endoscopic accessories and all connected devices with regard to functionality and any visual and mechanical damages in order to exclude risk of injury! Endoscopes with damages or defective individual parts may no longer be used. When using accessories and further components, the complete function and intended use of the medical products must be guaranteed and checked by the user before use.

If in doubt, contact your distributor or the manufacturer.

The device is not intended for operation in explosive areas!

Mechanical strain caused by dropping, over-bending, bending at tight angles, strong blows and torsional, tractional or pressure strains can lead to damage or the destruction of the endoscope and, thus, to impaired operability!

We do not accept liability for damages caused by wrong handling or misuse of the endoscope!

Endoscopic procedures should only be performed by trained specialists (e.g. physicians) with corresponding training in and familiarity with endoscopic procedures, whereby it is the responsibility of the user to continuously consult information regarding indications, contraindications, possible complications and risks as well as the development of the endoscopic procedure.

Attention: Incompatible and incorrectly combined products could injure the patient or put him at risk.

Information regarding the use with camera systems

Endoscopes are used to visualize the operating field. You need to be connected to a compatible camera and light source system for this.

joimax® combo endoscopes can only be used with joimax® Camera and Light Source systems (C-Camsource®, C-Camsource® HD Twister, etc.). The camera is connected to the endoscope by means of a cable via which light is produced and the optical image is transferred. The camera cable head and the endoscope are interconnected with a locking mechanism.

joimax® endoscopes with ocular cones can be adapted to all conventional camera and light source systems for medical endoscopes. Here, the camera cable is connected to the endoscope's ocular cone in accordance with ISO/TS 18339:2015 and the light cable is screwed onto the manufacturer-specific adapter.

Notes on use with light sources

The failure of a used light source may lead to hazards. Therefore you should have a replacement lamp readily available or use light sources with a replacement lamp. Do not move the endoscope while replacing the lamp. If possible, safely remove the endoscope.

Please note that light is an energy form that could heat up the endoscope tip.

When combined with high-performance light sources, both the optical fiber end of the light source as well as the instrument can reach temperatures that may result in burns. High radiation energy light can further lead to temperature increase in tissue. In the case of invasive application, it is imperative to avoid temperatures above +41°C, as this may lead to tissue damage. You should, thus, avoid direct tissue contact and, if possible, ensure sufficient rinsing of the operative field.

Do not stare into the light or onto surfaces reflecting the light beam. The intense light may cause temporary loss of vision.

Notes on use with flushing pumps

The joimax® endoscopes with working channel and the shafts of the joimax® nucleoscope have 1-2 flush connectors (Luer Lock) for flushing or aspiration. Only compatible tube sets may be connected to the flush connectors. joimax® tubing sets and patient lines (e.g. JTSB350D, JTSB300D and JPPL0350) are compatible with joimax® endoscopes. Only flushing solutions suitable for endoscopic procedures may be used for flushing. Always check whether the connection between endoscope and tubing is free of leakage before starting the procedure. More information is available in the joimax® surgery instructions. When using a suction pump: Do not to set the suction pressure higher than the irrigation pressure to avoid unintentional drop in pressure at the cavity.

Notes on use with sealing caps

Make sure the sealing cap has a tight fit and does not allow leakage of irrigation fluid other than through the intended opening.

Make sure the opening of the sealing cap is not blocked and irrigation fluid can flow through freely.

Notes on the combination with other medical products

The adaptation to other devices or accessories (such as shavers, insufflators, HF devices, working elements, laser devices, pneumatic or electro-hydraulic lithotripters and others) creates the most varied therapeutic application options, whereby it is imperative that you adhere to the device's and its accessories' instructions for use and safety instructions and that you ensure that the user is sufficiently trained.

If in doubt, contact your dealer or the manufacturer.

Protective measures for HF, laser and other high-energy applications are not integrated. Note that you may only adapt devices approved for medical purposes!

Detailed understanding of principles and methods applied in laser endoscopy and during electro-surgical procedures is required to prevent risk of shock and burns for patients and operators as well as damages of other devices and instruments.

We do not accept liability claims caused by incorrect use or the combination with other devices and instruments.

When operating an endoscope together with electro-medical devices, ensure the adherence to BF-conditions (isolated, ungrounded application part according to DIN ISO 60601-1).

The leakage currents may add up if endoscopes are used with electro-medical devices and/or energetically operated endoscopically usable accessories.

Notes on the use of high-frequency surgical devices

Before using endoscopic high-frequency surgical devices, patients should be suitably prepared for the intended operation. This especially includes measures to remove and prevent the creation of ignitable gases.

In contrast to conventional high-frequency surgery, unsuitable (especially too low) performance settings in endoscopic high-frequency surgery can cause significant penetration into the surrounding tissue.

Performance settings should, therefore, be made according to the operator's experience with reference to associated clinical references and/or suitable training. To prevent burns and/or undesired penetration into the surrounding tissue or damages to the endoscope, the high-frequency current should only be activated once the respective application part (electrode) can be seen through the endoscope. Never touch the endoscope with an activated electrode.

Notes on use with laser devices or laser probes with endoscopes

If endoscopes or endoscopic accessories are used together with laser devices, wear suitable protective eye-wear to prevent potential damage to the eyes. To prevent hazards and in relation to possible application limitations, follow the respective device-specific instructions for use and safety notes when performing with laser devices. Never touch the endoscope with an activated laser fiber.

Notes on use with lithotripters

To prevent hazards and in relation to possible application limitations, follow the respective device-specific instructions for use and safety notes when performing ultrasound, electro-hydraulic, pneumatic and mechanical lithotripsy.

CT (computer tomography)/X-ray

X-ray radiation may discolor optical components and, thus, damage endoscopes.

A simultaneous application of CT (computer tomography)/X-rays and endoscopes may lead to hazards.

MRI/Magnetic resonance imaging

Endoscopes containing stainless steel are not suitable for MRI.

Adverse Effects

The following list shows a non-exhaustive overview of possible adverse effects that may be associated with endoscopic spinal surgery:

- Reherniation
- Transient dysesthesia
- Headaches due to dural tear
- Transitory foot drop
- Transitory sensibility disturbance of the foot
- Transient numbness
- Cauda equina syndrome
- Dural tear
- Infection
- Thrombophlebitis
- Vascular injury

4. SCOPE OF INSTRUCTIONS

This document applies to the following joimax® products:

4.1 ENDOSCOPES WITH WORKING CHANNEL

REF	Description	Compatible with #
FS5831120C	TESSYS® Foraminoscope thoracic, combo	4
FS58311200	TESSYS® Foraminoscope thoracic ocular	
FS6342181C	TESSYS® Foraminoscope, combo	1
FS63421810	TESSYS® Foraminoscope, ocular	
FS7347171C	TESSYS® Foraminoscope, combo	3
FS73471710	TESSYS® Foraminoscope, ocular	
FX6342208C	TESSYS® Foraminoscope XT, combo	2
FX63422080	TESSYS® Foraminoscope XT, ocular	

REF	Description	Compatible with #
LS1006125C	iLESSYS® Delta Laminoscope, combo	6
LS1006125O	iLESSYS® Delta Laminoscope, ocular	
LS6342125C	iLESSYS® Laminoscope, combo	7
LS6342125O	iLESSYS® Laminoscope, ocular	
LS7347068C	CESSYS® Dorsal Laminoscope, combo	9
LS7347068O	CESSYS® Dorsal Laminoscope, ocular	
LS7347150C	iLESSYS® Pro Laminoscope, combo	5
LS7347150O	iLESSYS® Pro Laminoscope, ocular	
MS5830125C	Multiscope, combo	8
MS5830125O	Multiscope, ocular	
HYBRID SCOPES		
CH402100C	CESSYS® Cervical Hybrid Scope, combo	10

COMPATIBLE PRODUCTS

#	REF	Instrument Set
1	TESBASIC	TESSYS® Transforaminal Endoscopic Surgical System
2	TESBASICXT	TESSYS® XT Transforaminal Endoscopic Surgical
	TESXTD	TESSYS® XT Extension Set for TESBASIC
3	TESSAP	TESSYS® Transforaminal Endoscopic Surgical System
4	TESTHX	TESSYS® Thorax Endoscopic Surgical System
5	ILESSYSPRO	iLESSYS® Pro Interlaminar Endoscopic Surgical System
6	ILESSYSDL	iLESSYS® Interlaminar Endoscopic Surgical System
7	ILESSYSSH	iLESSYS® Interlaminar Endoscopic Surgical System
	MUSFV	MultiZYTE® Endoscopic Surgical System
	MUSBASIC	MultiZYTE® Endoscopic Surgical System
8	MUSCER	MultiZYTE® Facet Cervical Surgical System
9	CESSYSDORS	CESSYS® Dorsal Endoscopic Surgical System
10	CESSYSANT01	CESSYS® Instrument Set
	CESSYSANT02	CESSYS® Instrument Set, oval

4.2 ENDOSCOPES WITHOUT WORKING CHANNEL

REF	Description	Compatible with #
PS1909132C	Nucleoscope, combo	12
PS1909132O	Nucleoscope, ocular	
PS2730246C	Nucleoscope, combo	11
PS2730246O	Nucleoscope, ocular	
BS4000175O-4K	Biportal Laminoscope, ocular	13
BS4015175O-4K	Biportal Laminoscope, ocular	
BS4030175O-4K	Biportal Laminoscope, ocular	

COMPATIBLE PRODUCTS

#	REF	Instrument Set
11	INTENTS01	intENTS® interv. Endoscopic Nucleus Therapy Set
12	INTENTSZ01	intENTS® interv. Endoscopic Nucleus Therapy Set
13	BSHS04175	Endoscope Shaft

5. INTENDED USE

5.1 INTENDED USE/INTENDED PURPOSE

joimax® endoscopes are optical devices enabling the observation of the body's interior through artificial cavities during endoscopic spinal surgery.

Additionally, joimax® endoscopes with working channel serve as guidance for other medical devices (e.g. surgical instruments) to the surgical site for the purpose of manipulation (e.g. palpation, removal, ablation) of tissue excluding all soft tissue structures of the central nervous system, namely spinal cord and meninges during minimally invasive spinal surgery. joimax® endoscopes with working channel are equipped with irrigation connectors. The irrigation connectors are intended to irrigate the surgical site with an irrigation pump.

Descriptions: The name of joimax® endoscopes (e.g. Foraminoscope®, Laminoscope®) are used to indicate the affiliation to an access approach. The intended use is identical.

5.2 INDICATIONS AND CONTRAINDICATIONS

LUMBAR AND THORACIC DEHERNIATION/DECOMPRESSION VIA TRANSFORAMINAL ACCESS (TESSYS® METHOD)

Indications

Indications for endoscopic spinal surgery in case of persisting clinical symptoms that could not be improved by conservative therapy:

- Herniated intervertebral disc (protrusion, prolapse, sequester and free fragment)
- Spinal stenosis through soft and/or bone tissue
- Degenerative disease of the ligamentum flavum
- Cauda equina syndrome
- Facet joint cyst

Contraindications for Lumbar Application

- L5/S1 with high iliac crest blocking the transforaminal access

Contraindications for Thoracic Application

- Heavily calcified disc herniation
- Ossification of the posterior longitudinal ligament (OPLL)
- Highly migrated disc herniations
- Severe disc narrowing

LUMBAR AND CERVICAL DEHERNIATION/DECOMPRESSION VIA INTERLAMINAR ACCESS (iLESSYS® METHOD)

Indications

Indications for endoscopic spinal surgery in case of persisting clinical symptoms that could not be improved by conservative therapy:

- Herniated intervertebral disc (protrusion, prolapse, sequester and free fragment)
- Spinal stenosis through soft and/or bone tissue
- Degenerative disease of the ligamentum flavum
- Cauda equina syndrome
- Facet joint cyst

Contraindications

- Major posterior scarring and adhesion resulting in unclear imaging of the anatomical structures at the surgical site

CERVICAL DEHERNIATION/DECOMPRESSION VIA ANTERIOR INTRADISCAL ACCESS (CESSYS® VENTRAL METHOD)

Indications

Indications for endoscopic spinal surgery in case of persisting clinical symptoms that could not be improved by conservative therapy:

- Herniated intervertebral disc (protrusion und prolapse, sequester)
- Mild spinal stenosis through soft and/or bone tissue

Contraindications

- Disc height smaller than 4 mm
- Calcified disc herniations
- Instability

CERVICAL DEHERNIATION/DECOMPRESSION VIA INTERLAMINAR ACCESS (CESSYS® DORSAL METHOD)

Indications

Indications for interventional endoscopic spinal surgery in case of persisting clinical symptoms that could not be improved by conservative therapy:

- Herniated intervertebral disc (protrusion, prolapse, sequester and free fragment)
- Spinal stenosis through soft and/or bone tissue
- Facet joint cyst

CERVICAL AND LUMBAR INTRADISCAL PAIN TREATMENT VIA DORSOLATERAL ACCESS (intENTS® METHOD)

Indications

Indications for interventional endoscopic spinal surgery in case of persisting clinical symptoms that could not be improved by conservative therapy:

- Discogenic pain
- Annular fissures (high intensity zones)
- Herniated intervertebral disc (protrusion and contained herniated disc)

Contraindications

- Sequestered herniation and free fragments
- Bony Stenosis
- Instability

LUMBAR AND CERVICAL FACET JOINT TREATMENT AND SIJ-TREATMENT (MultiZYTE® METHODS)

Indications for endoscopic spinal surgery in case of persisting clinical symptoms that could not be improved by conservative therapy:

Indications Lumbar/Sacroiliac Joint (SIJ)

- Chronic Low Back Pain
- Degenerative disease of the facet joint
- Facet joint hypertrophy
- Inflammation of the facet joint and/or SIJ

Indications Cervical

- Cervicogenic Headache
- Chronic neck pain
- Degenerative disease of the facet joint
- Inflammation of the facet joint

Contraindications

- Instability

GENERAL CONTRAINDICATIONS (ALL METHODS)

The application of joimax® instruments is contraindicated if minimal invasive spinal surgery is contraindicated. Contraindications for diagnostic and endoscopic spinal surgery include:

- General surgical contraindications
- Physiologically or psychologically inadequate patient
- Infections
- Patients with inadequate blood coagulation

Based on the general condition of the patient, the responsible physician/surgeon must decide by interdisciplinary consideration whether the planned operation is feasible or not. The laws and regulations applicable in the country of use must be complied with. These contraindications can be relative or absolute. The surgeon must assess the benefits to the risks and inform the patient about his or her decision. The above lists are not exhaustive.

5.3 INTENDED PATIENT POPULATION

The instruments are intended to be used on grown up patients of all genders with appropriate indications.

5.4 INTENDED USERS/USER POPULATION

The instruments are intended to be used only by certified medical professionals who have been specifically trained (e.g. orthopedic surgeons, neurosurgeons) and who are familiar with respective surgical techniques.

5.5 AREAS OF APPLICATION

The instruments are intended to be used in a designated clinical facility. The instruments are handled by a physician in a sterile environment.

6. CONDITION OF DELIVERY

Carefully remove the endoscope as well as its accessories from the packaging and remove all packaging components.

Delivery contains:	1	Endoscope
	1	Instructions for use
	2	Cleaning brush for endoscope type with irrigation and working channel
	1	Cleaning adapter for Foraminoscope®, Laminoscope®, Multiscope

Immediately after receipt check the endoscope and its accessories for completeness and obvious damages. Transport damages can only be claimed if the supplier is informed immediately (within 24 hours) after delivery.

Always use the original packaging for any return consignment of devices or accessories. Describe the faults or malfunction and state a contact person for any enquiries.

Ensure that only sterile endoscopes are returned.


Caution: When connecting accessories and other medical equipment, please follow the notes on combination with other medical products listed in chapter 4. "Scope of instruction".


The endoscopes are delivered in unsterile condition and must be cleaned and disinfected or sterilized before their initial use in accordance with the reprocessing instructions (see chapter 9.1 "Cleaning and Disinfection" and chapter 10, "Sterilization").

7. PRODUCT DESCRIPTION (ASSEMBLY / SKETCH / DESCRIPTION OF PARTS)


Example Nucleoscope (Double-cable-technique (Ocular))

Nucleoscope Shaft


Nucleoscope Shaft-Obturator


Example Arthroscope


(Single-cable-technique (Combo))

Example Foraminoscope® / Laminoscope® / Multiscope (Single-cable-technique (Combo))


Example Foraminoscope® / Laminoscope® / Multiscope (Double-cable-technique (Ocular))


Example cervical hybrid endoscope

(Single-cable-technique (Combo))

Cleaning Adapter for Foraminoscope® / Laminoscope® / Multiscope

Just mount the cleaning adapter onto the end of the proximal working channel.

Cleaning Adapter with Luer Lock

Compatibility

- Optical fiber connection
 - Imaging device connection
- Compatible with prevalent systems
Ocular cone, screwable according to ISO/TS 18339:2015

Biocompatibility

All materials coming into body contact are biocompatible.

8. OPERATING OF THE DEVICE

The installation of the minimum equipment required for endoscopy (this can vary through the combination of different devices):




The following pictures refer to the handling of the equipment under non-sterile conditions.

Example Foraminoscope® / Laminoscope® / Multiscope

a) Single-cable-technique (Combo)

- 1) Align light and optical channel of the endoscope and camera head to each other
- 2) Pull back bayonet lock (small lever)
- 3) Insert the endoscope to stop in the camera head coupling
- 4) Push bayonet lock (small lever) back in starting position to ensure locking of the coupling
- 5) Turn small lever to the right to fixate endoscope
- 6) Manually focus the endoscopic image using the big lever


Please note the operating instructions of the connected camera.

b) Double-cable-technique (Ocular):

TV adapter connection to the endoscope e.g. to an arthroscope


Ensure that the endoscope connection is locked

Connection of fiberoptical cable



Ensure that the fiber optic cable is locked

If applicable: Instrument connection

Ensure locking of the instrument

Please note the operating instructions of the connected camera.

9. CLEANING AND DISINFECTION

Preliminary remarks

The following suggested procedure with regard to preparation and care must be regarded as recommendation. Endoscopes can also be disinfected and cleaned as well as sterilized in other ways.

The operator is responsible for any preparation procedures with regard to achieving the desired disinfection, cleaning and sterilization effect. Also keep in mind that the type and manner of preparation can have an influence on the durability of the endoscope.

Check your endoscope for functionality (such as sufficient lighting of fibers, a clear, focused, bright and round image) and for possible damages (such as sharp edges, visible material deformation or loose parts) before each use.

The endoscope has to be cleaned and sterilized before the first use. The endoscope should be cleaned immediately after each application. Preferably use demineralized water for cleaning in order to prevent damaging effects to the components of the endoscope. Sterilize the endoscope before each further use or disinfect in accordance with the medical indication.

 **Caution: Inadequate or omitted reprocessing can lead to infection or inflammation of the patient.**

The preparation procedure recommended below conforms with the DIN EN ISO 17664:2004 on the sterilization of medical devices – the preparation of resterilizable medical devices. The procedures are validated. Irrespective of the above-mentioned, the user is responsible for achieving a sterilization effect.

9.1 REPROCESSING INSTRUCTIONS

Rinsing and suction valves, when present, should always be dismantled before sterilization. If possible, use a automated procedure (disinfector) for cleaning and disinfection. Due to uncertain reproducibility, manual procedures should only be used if a automated procedure is unavailable.

9.2 RESTRICTIONS AND LIMITATIONS IN REPROCESSING, DURABILITY

Frequent reprocessing of endoscopes affects their functionality. The end of the product life cycle is usually determined by wear and damages through use.

Should an inspection (refer to chapter 9.7. "Control, maintenance and inspection") reveal defects, return the endoscope to the manufacturer for evaluation or repair.

The reprocessing of the endoscopes has been proven for up to 100 cycles. Depending on individual handling, application and reprocessing, the actual figure may differ.

9.3 HANDLING AFTER USE

Endoscopes must be detached from other devices and cleaned by wiping immediately after use. Remove surface dirt with a lint-free, soft cloth / paper towel. Should instruments and endoscopes have come into contact with corrosive etching and medicinal agents during application, they must be cleaned immediately.

"Throwing" instruments (endoscopes) into a collecting container during an operation will invariably lead to damages. Ensure that instruments are deposited gently after use.

9.4 PREPARATION BEFORE CLEANING

In order to prevent blood, proteins or other substances drying on the endoscope and to protect the staff, dispose of the endoscope in wet condition immediately after surgical use. Dried proteins on the endoscope complicate cleaning, disinfection and sterilization.

Caution: Endoscopes may not be decontaminated and cleaned in an ultrasonic bath.

Caution: Use approved cleaning agents only. Do not submerge the endoscope in alcohol or other aggressive liquids. Hydrogen peroxide solution (H_2O_2 solution) is not tested yet and, therefore, not fully recommendable.

Procedure

Use running water or a disinfectant solution; the disinfectant should be free of aldehydes (these fixate blood residues), have a certified effect (e.g. be approved by VAH/DGHM or FDA or bear a CE mark), be suitable for the disinfection of instruments and be compatible with the instruments. To manually remove dirt, only use a soft brush or a clean, soft cloth used for this purpose only.

Never use metal brushes, steel wool or other sharp objects.

Rinse all instrument lumina five times using a disposable syringe (minimum volume 10-50 ml (in accordance with the instrument size)). Disassemble the rinsing valves (does not apply to endoscopes without rinsing valves). Disconnect the endoscope from the connections (rinsing and suction tube, light conducting cable, camera head). If equipped, remove any telescope bridges or stop-cock valves from the endoscope (does not apply to endoscopes without rinsing valves).

Remove any light conducting adapters and screw caps on valves, stop-cock keys and protective covers on valves and bridges. Remove coarse surface dirt with a moist, lint-free, soft cloth.

Remove all sealing caps, valves and duct nozzles from the endoscope.

9.5 CLEANING

Cleaning is required for hygienic reasons.

9.5.1 AUTOMATED CLEANING / DISINFECTION (DISINFECTOR / CDD (CLEANING AND DISINFECTING DEVICE))

Endoscopes with working channels can be cleaned in respective washers with special endoscopic cleaning programs, but the working and rinsing channels must first be manually cleaned as described in the next chapter, "manual cleaning". Connect the working and rinsing channels to the machine's rinsing connection. To connect the working channel to the machine's rinsing connection, use the supplied rinsing adapter depending on the endoscope version.

⚠ In the case of thermal cleaning procedures, ensure that cooling periods are adhered to. Accelerating the cooling process may lead to damaged optical components!

The respective automated processes must ensure that the endoscopes remain safely positioned on the instrument carriers and do not detach or cause damage to one another.

When selecting the disinfectant, make sure that:

- the disinfectant has proven effectiveness (e.g. DGHM or FDA approval or CE marking according to DIN EN ISO 15883)
- if possible, a tested program for thermal disinfection (AO value > 3000 or - for older devices - at least 5 min at 90° C/194° F) is used
- the used program is suitable for the instruments (see chapter 10.6 "Reusability") and contains sufficient rinsing cycles
- for rinsing only sterile or low-germ (max 10 germs/ml) and endotoxin-poor (maximum 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water) is used
- the air used for drying is filtered
- the disinfectant is regularly maintained and checked

When selecting the applied cleaning system, ensure

- that it is suitable for cleaning metal and plastic instruments,
- that – unless thermal disinfection is used – a suitable disinfectant with proved effectiveness (e.g. VAH /DGHM or FDA approved or with a CE mark) is used and that this is compatible with the applied cleanser,
- that the applied chemicals are compatible with the instruments (refer to chapter 10.5 "Material durability").

It is imperative to adhere to concentrations specified by the manufacturer of the cleanser or disinfectant.

Procedure:

1. Disconnect the endoscope to the greatest possible extent (light conductor connection, handle, possibly rinsing valves).

2. **Place the disassembled instruments into the disinfectant, ensuring that the instruments do not touch one another. If applicable: Connect all instrument lumina to the rinsing connection of the disinfectant using the available LuerLock connection or a suitable rinsing adapter.**

3. Start the program.
4. After the program has completed, remove the instruments from the disinfectant.
5. Inspect and pack the instruments, preferably immediately after removal (refer to chapter 9.7 "Control, maintenance and inspection" and chapter 10.1 "Packaging") or after additional drying in a clean place, if applicable.

Proof of fundamental instrument suitability for effective automated cleaning and disinfection was provided by an independent accredited testing laboratory with the use of the disinfectant G 7836 CD (thermal disinfection, Miele and Cie. GmbH & Co., Gütersloh) and the cleaning agent neodisher® MediClean forte (Dr. Weigert GmbH & Co. KG, Hamburg), whereby the procedure described above was used.

9.5.2 MANUAL CLEANING AND DISINFECTION

When selecting the applied cleansers and disinfectants, ensure

- that they are suitable for cleaning metal and plastic instruments,
- that a suitable disinfectant with tested effectiveness (e.g. VAH / DGHM or FDA approved or with a CE mark) is used and that this is compatible with the applied cleanser.
- that the applied chemicals are compatible with the instruments (refer to chapter 10.5 "Material durability").

If possible, avoid combined cleansers / disinfectants. Combined cleansers / disinfectants should only be used in cases with very little contamination (no visible dirt).

It is imperative to adhere to concentrations and reaction times specified by the manufacturer of the cleanser or disinfectant. Only use recently made solutions, sterile or germ-reduced (max. 10 germs / ml) as well as endotoxin reduced (max. 0.25 endotoxin units / ml) water (e.g. purified / highly purified water) or filtered air for drying.

Procedure:

Cleaning

1. Disconnect the endoscope to the greatest possible extent (light conductor connection, handle, possibly rinsing valves).
2. Submerge the disassembled instruments into the cleaning bath for the specified reaction time so that the instruments are sufficiently covered (careful brushing with a soft brush, if required), ensuring that the instruments do not touch one another. Move movable parts around several times during cleaning. If applicable: At the beginning and at the end of the reaction time, rinse all instrument lumina five times by using a disposable syringe (minimum volume 10-50 ml (in accordance with the instrument's size)) and the LuerLock connection or a suitable rinsing adapter.
3. Then remove the instruments from the cleaning bath and thoroughly rinse them with water at least three times. If applicable: At the beginning and at the end of the reaction time, rinse all instrument lumina five times by using a disposable syringe (minimum volume 10-50 ml (in accordance with the instrument's size)) and the LuerLock connection or a suitable rinsing adapter.
4. Check the instruments (refer to chapter 9.7 "Control, maintenance and inspection").

Disinfection

5. Submerge the disassembled, cleaned and checked instruments into the disinfection bath for the specified reaction time so that the instruments are sufficiently covered, ensuring that the instruments do not touch one another. Move moveable parts around several times during disinfection. If applicable: At the beginning and at the end of the reaction time, rinse all instrument lumina five times by using a disposable syringe (minimum volume 10-50 ml (in accordance with the instrument's size)) and the LuerLock connection or a suitable rinsing adapter.
6. Then remove the instruments from the cleaning bath and thoroughly rinse them with water at least five times. If applicable: At the beginning and at the end of the reaction time, rinse all instrument lumina five times by using a disposable syringe (minimum volume 10-50 ml (in accordance with the instrument's size)) and the LuerLock connection or a suitable rinsing adapter.
7. Dry the instruments by blowing off and out with filtered pressurized air.
8. Pack the instruments, preferably immediately after removal (refer to chapter 10.1 "Packaging", after additional drying in a clean place, if applicable).

Proof of fundamental instrument suitability for effective manual cleaning and disinfection was provided by an independent accredited testing laboratory with the use of the cleanser Cidezyme/Enzol and the disinfectant Cidex® OPA (Johnson & Johnson GmbH, Norderstedt), whereby the procedure described above was used.

9.6 DRYING – GENERAL REQUIREMENT

Endoscopes and their working and rinsing channels or cavities must be completely dry and then stored in suitable containers.

9.7 CONTROL, MAINTENANCE AND INSPECTION

After cleaning or cleaning and disinfection, check all instruments for corrosion, chipping and soiling and discard damaged instruments (numerical restrictions of reuse refer to chapter 10.6 "Reusability").

Even if carefully cleaned and disinfected, optics may still have adhesions on their mechanisms and surfaces. Use pure alcohol to clean the optics glass surfaces in order to ensure perfect image quality.

After the cleaning it is recommended to grease the rinsing valves (not applicable for endoscopes without rinsing valves). For this, use only special grease for rinsing valves. An ongoing maintenance of the rinsing valves is not required.

Checking the mechanism and endoscope surface

The endoscope surfaces must be undamaged and especially free of any sharp edges. Note any dents, kinks, mechanical/thermal damages, e.g. caused by HF or laser-operated surgical devices as well as any tears or bursts. Damages of the endoscope can lead to loss of function. Damages such as sharp edges or protruding burrs can injure the patient.

Checking the fiber optics

Hold the distal end of the endoscope towards a lit window or a bright light or similar.

Caution: Do not use a cold-light source for this test! Looking directly into the radiated light of a cold light source can lead to eye damage.

Examine the optical fiber connection. The individual fibers should appear brightly. Point the distal end of the scope towards a bright light and move it slightly up and down. The brightness of the fibers at the proximal end will change slightly. There is no need to worry if individual fibers remain dark, but if more than 20 to 30% of fibers are broken, then endoscopic work is extremely handicapped.

The surfaces of the light entrance and light exit areas should be smooth and clean. Deposits or individual fibers, present on the surfaces, which are coarse and/or palpable, can lead to insufficient illumination. If the surfaces are not cleaned, then further use or preparation can lead to progressive damage to the endoscope.

Caution: Endoscopes with damaged fiber optics should be returned to the manufacturer or an authorized service agent for evaluation and, if applicable, for repair.

Checking the proximal and distal glass surfaces

- The glass surfaces must be clean and free of deposits. If the visual check of the cover glass reveals tenacious or dried contaminations, these can be removed with suitable cleaning pastes or with the help of an alcohol soaked cotton swab or toothpick. A frequent cause for these deposits is insufficient rinsing of the optics after cleaning and disinfection.
- The image must be focused and clear. An unfocussed, unround, blurred, foggy image is a sign of damage.

Caution: If cleaning does not remove persistent deposits, the endoscope should be sent to the manufacturer or an authorized service agent for inspection.

Caution: Endoscopes with damaged glass surfaces (e.g. chipping), impaired image quality or conspicuous surface damage and deformation may not be used further and should be discarded or returned to the manufacturer or an authorized service agent for inspection.

Maintenance

Reassemble the disassembled instruments (light conductor connection, handles, possibly rinsing valves).

Do not apply instrument oils. Lubricating the rinsing valves is only permitted after sterilization and with medical grease with proven biocompatibility (does not apply to endoscopes without rinsing valves).

Regular maintenance by the manufacturer is not required.

10. STERILIZATION

 **Caution:** It is not advisable to sterilize insufficiently cleaned endoscopes as the success of sterilization is extremely dependent on a thorough pre-cleaning! Before sterilization, make sure the endoscope is clean – especially the optical surfaces – and check that the endoscope is – without restrictions – ready for use (see chapter 9.7 "Control, maintenance and inspection").

 **The user is responsible for performing the specified sterilization process with regard to achieving the desired or required sterilization effect.**

 **Caution:** Inadequate or omitted reprocessing can lead to infection or inflammation of the patient.

Clean any optical glass with pure alcohol prior to sterilization!

During sterilization, ensure that optical components do not come into contact with hot metal areas, as these thermal bridges can lead to the destruction of the applied materials and, thus, to a leakage of the entire system.

10.1 PACKAGING

The packaging of the sterile goods provides protection of the content during transportation. The material must conform to the required characteristics. The durability against steam, gas, plasma as well as temperature differences and pressure changes must be stated.

Sort the cleaned and disinfected instruments into the provided sterilization tray.

Please pack the instruments or the sterilization trays into disposable packaging (disposable or dual packaging) and/or sterilization containers that comply with the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607-1/11607-2
- suitable for steam sterilization (heat resistant minimum 141 °C (286 °F) sufficient vapor permeability)
- sufficient protection of the instruments or sterilization packaging from mechanical damage
- regular maintenance according to manufacturer specifications (sterilization container)

Packaging must be marked with the sterilization type and date, expiry date and batch number.

 **Caution:** Sterilized endoscopes must always be stored in sterile goods packaging or storage containers within a closed area or designated cabinet, safely protected from heat, radiation, dust, moisture and temperature fluctuations.

10.2 STEAM STERILIZATION/AUTOCCLAVING

joimax® endoscopes are suitable for sterilization with a temperature of 132/134 °C max. with a holding time of 5 minutes (max. 18 min.). This creates an overall pressure of 3.1 bar in the autoclave, which corresponds to an overpressure of approx. 2.2 bar on the autoclave display. This value of 2.2 bar is indicated also on the endoscopes.

max. 134°C
 AUTOCLAVE
 max. 2.2 bar / 18 min.

Marking on the endoscopes

Information about validation and ISO standard

- validated in accordance with DIN EN ISO 17665-1:2006 – (valid IQ / OQ (commissioning) and product-specific provision of services (PQ))
- maximum sterilization temperature 138 °C (280 °F; plus tolerance in accordance with DIN EN ISO 17665-1:2006)
- steam sterilizer in accordance with DIN EN 13060 or DIN EN 285

Sterilization method	Temp. °C	Temp. °F	Temp. holding time	Drying time
Pre-vacuum	121 °C	250 °F	20 min.	20-30 min.
	132 °C	270 °F	3 min.	20-30 min.
	134 °C	273 °F	3 min.	20-30 min.
	135 °C	275 °F	3 min.	20-30 min.
Prion program	134 °C	273 °F	18 min.	20-30 min.

Attention:

- The high-speed sterilization procedure is not permitted.
- Never use hot-air-sterilization, radiation sterilization, formaldehyde or ethylene-oxide sterilization procedures.
- Other autoclavation settings and cycles can have a negative impact on the device or its components.
- Do not exceed the weight and loading of the goods to be sterilized as this may lead to excess condensate and, thus, to damages caused by rust.
- When sterilizing several instruments in one sterilization cycle, do not exceed the maximum load of the sterilizer.
- Please note that sterilization success is decreased with each increase in the number of instruments in one sterilization cycle. Please read the manufacturer specifications.
- The endoscope must be dried after cooling to room temperature.
- Adhere to the specifications of the sterilizer manufacturer, especially with regard to aeration times after sterilization. Observe the valid national legal regulations.

Always adhere to cooling periods specified for thermal sterilization procedures. Accelerating the cooling process may lead to damaged optical components!

Strictly follow the instructions for use of the sterilizer's manufacturer.

10.3 PLASMA STERILIZATION

The endoscopes may be sterilized with the following STERRAD® System manufactured by ASP, a Johnson & Johnson company:

- STERRAD® 100S

Caution: It is mandatory to follow the instructions for use provided by the manufacturer of the STERRAD® System.

The effectiveness of the STERRAD® 100S System by using the “Short cycle with booster” has been proven.

For the plasma sterilization procedures, we recommend a periodical maintenance of the scopes in order to prolongate the life cycle of the scopes.

For maintenance or repair request, please return the scopes to the manufacturer or your local distributor.

10.4 STORAGE

The endoscope must have been completely dried.

Always store sterilized endoscopes in packaging or storage containers suitable for sterile goods in a closed sterile area or designated cabinet, protected from heat, radiation, humidity and temperature fluctuations as well as contamination.

10.5 MATERIAL DURABILITY

When selecting cleansers and disinfectants, always ensure that they do not contain the following components:

- Organic, mineral or oxidating acids (minimum admissible pH-value 5.5)
- Strong lyes (maximum admissible pH-value 11; we recommend a neutral/enzymatic or slightly alkaline cleanser)
- Organic solvents (e.g. alcohol, ether, ketones, benzines)
- Oxidants (e.g. hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons
- Oils

Never clean instruments, sterilization trays and sterilization containers with metal brushes or steel wool.

All instruments may only be exposed to temperatures below 138 °C (280 °F; plus tolerance in accordance with DIN EN ISO 17665-1:2006)!

10.6 REUSABILITY

Endoscopes can – if carefully handled and if undamaged and unsoiled – be reused multiple times. Every subsequent use or the use of damaged and/or soiled instruments is the responsibility of the user.

We refuse any liability in the event of non-observance.

Steam sterilization/Autoclaving:

Autoclaving at 134 °C (273 °F) and holding the temperature for 18 minutes (Prion program) results in a faster aging of the endoscope as compared to holding the temperature for 5 minutes only. joimax® recommends to return the scope periodically to the manufacturer for precautionary maintenance if 18 minutes autoclaving is used.

ASP STERRAD®:

joimax® recommends to return the scope periodically to the manufacturer for precautionary maintenance when using the STERRAD® 100S.

10.7 ADDITIONAL INFORMATION

The manufacturers have validated the instructions listed above as suitable for medical devices serving the preparation of a medical devices for the purpose of its reuse. The processor is responsible for ensuring that the actually performed preparation achieves the desired results with the used equipment, materials and staff at the processing facility. This normally requires validation and routine monitoring of the process. Any deviations from the provided instructions performed by the processor must be carefully assessed with regard to their effectiveness and possible disadvantageous consequences.

11. TROUBLESHOOTING

Problem	Possible cause	Remedy
Image is blurred, foggy	Dirty glass surfaces	Clean glass surfaces acc. to chapter 9.5.2 (manual cleaning and disinfection)
	Tenacious or dried coating on the glass surfaces	Remove coating acc. to instructions under chapter 9.5.2 (manual cleaning and disinfection), Check water quality
	Leaking, defective lens system	Send endoscope in for repair
Image too dark, insufficient illumination	Dirty glass surfaces	Clean glass surfaces acc. to chapter 9.5.2 (manual cleaning and disinfection)
	Tenacious or dried coating on glass surfaces	Remove coating acc. to instructions under chapter 9.5.2 (manual cleaning and disinfection), check water quality
	Wrong light conducting cable connection	Check seating of light conductor cable, exchange if necessary
	Defective fiber optics	Check fiber optic acc. to chapter 9.7 (control, maintance and inspection)
	Defective light conductor cable(s) Light source	Check light conductor cable, light source

Yellow-tinted illumination	Dirty fiber optics	Clean glass surfaces acc. to chapter 9.5.2 (manual cleaning and disinfection), send endoscope to service, if required
	Dirty, defective shine light conducting cable	Check light conducting cable (e.g. onto a white surface), replace if necessary
Spot formation, discoloration	Insufficient cleaning (e.g. residual proteins)	Clean again, if necessary, with thorough wiping
	Insufficient rinsing of the endoscope between phases (especially prior to sterilization)	Ensure sufficient rinsing preparation (see chapter 9.5.1. (automated cleaning/disinfection) and chapter 9.5.2 (manual cleaning and disinfection))
	Contaminated, excessive use of disinfectants and cleansers	Renew regularly the cleaning solutions
Leaky	Connections	Check connections between sealing cap, possibly rinsing valves
	Defective rinsing valves	Send endoscope in for repair

12. DISPOSAL

The disposal of medical devices and packaging must comply with the laws and regulations for waste handling in the respective country. Contact your local waste management company for further information.

13. SERVICE AND REPAIRS

In case of required maintenance or repair, please contact the joimax® GmbH service department or an authorized service agent.

In the interest of rapid processing of service enquiries, please send us the product with the following information:

- Article number (REF)
- Serial number (SN)
- Description of the error as accurate as possible.

Caution: Clean and sterilize the endoscope (incl. accessories) thoroughly before dispatch to protect your personnel as well as joimax® staff members. Should this be impossible for urgency reasons, prepare the endoscope to the greatest possible extent and mark accordingly.

For safety reasons, joimax® GmbH can refuse soiled or contaminated products.

All guarantee and warranty claims are lost if the user himself or an unauthorized repair company performs maintenance and repairs.

Arbitrary opening, repairing or modifying the device release the manufacturer of any liability for the operating safety of the system. No liability is assumed during the guarantee period as a result of such action.

The manufacturer is obligated to a warranty period of 12 months with regard to the device’s function. This warranty is limited to claims received within the warranty period, starting with the purchase of the device. This warranty only applies to defects which cannot be attributed to normal wear and tear, misuse or incorrect handling, lack of care or force majeure. This warranty excludes wear parts.

Always specify type and serial number when making enquiries or ordering spare parts.

We will be glad to answer any questions you might have.

English

13.1 EMERGENCY NUMBER

Please notify us of any events according to MPG § 29 or § 2 of the Medical Devices Safety Plan Ordinance under the following telephone number:

for the US 949-859-3472 (24-hour availability)

for all other countries +49 (0) 721 255 14-999 (24-hour availability).

Thank you for your support!

14. SIGNS AND SYMBOLS

	Please note!
	Caution/Attention
	Consult instructions for use
	Follow instructions for use
	CE marking in compliance with Directive MDD 93/42/EEC (until Dec. 23th 2023) and with the Medical Device Regulation 2017/745 (from Dec. 24th 2023 onwards)
	Serial Number
	Catalogue number
	Diameter
	Length
	Manufacturer
	Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician

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joimax®, Inc.