joimax®
Endoscopic Spine Experts

English

joimax® Vaporflex® Handle – Instructions For Use

Deutsch

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joimax® Vaporflex® 手柄 – 操作说明

TABLE OF CONTENT 2

GENERAL INFORMATION	3
PRODUCT DESCRIPTION	3
INDICATIONS OF USE	3
CONTRAINDICATIONS	3
TECHNICAL DATA	4
COMPONENTS	4
ASSEMBLY INSTRUCTIONS	4
DISASSEMBLY INSTRUCTIONS	5
FUNCTION TEST	5
MATERIAL DURABILITY	5
PRODUCT LIFETIME	5
NOTES ON THE CORRECT USE OF THE DEVICE	6
WARNINGS AND POTENTIAL RISKS	6
CAUTIONS	6
ELECTRO MAGNETIC COMPATIBILITY (EMC)	7
PRINCIPLE OF BIPOLAR ELECTROSURGERY	7
LIABILITY AND WARRANTY	7
LIMITED WARRANTY	7
PRIOR TO USE	8
INSTRUCTIONS FOR USE	8
PROCESSING INSTRUCTIONS	8
PRETREATMENT	8
MANUAL CLEANING AND DISINFECTION – FDA CLEARED	9
CLEANING	9
DISINFECTION	9
AUTOMATED CLEANING AND DISINFECTION – NOT FDA CLEARED, FOR OUTSIDE US MARKETS ONLY	10
PACKING	10
STERILIZATION	11
STORAGE	11
CUSTOMER SERVICE	11
SIGNS AND SYMBOLS	12
ORDER INFORMATION	12

GENERAL INFORMATION

The following instruction is only to demonstrate the correct set-up of the joimax® Vaporflex® Handle and does not give any recommendation for the surgical use of the instrument.

PRODUCT DESCRIPTION

The joimax® Vaporflex® Handle is an electrosurgical instrument for use in minimally invasive procedures. joimax® offers different bipolar probes varying in length, diameter and tip configuration. Various shaft lengths and/or shaft versions are available depending on the procedure and length of the probe.

The probes are available from joimax® as individual sterile packaging. The Vaporflex® Handle comprises different components, available either in a set or individually. These are non-sterile and must be reprocessed and sterilized before the first use according to the reprocessing instructions.

The Vaporflex® instrument comprises: handle, shaft and RF/HF-cable. The device must be fitted with a suitable disposable bipolar probe.

The joimax® Vaporflex® instrument is designed and intended for use only with electrosurgical RF/HF generators that have been tested at least to IEC 60601-1 3rd ed., IEC 60601-2-2, IEC 60601-1-2 (i.e. elliquence® SurgiMax®, SurgiMax® Plus, SurgiMax® Vapor, joimax® Endovapor®², joimax® Endovapor®)¹. Refer to the RF/HF generators documentation to ensure compatibility.

The joimax® Vaporflex® Handle is a Class I medical device in accordance with MDD 93/42/EEC.

INDICATIONS OF USE

The joimax® Vaporflex® and Legato® bipolar electrosurgical instruments are used for coagulation of soft tissue during open or minimal invasive surgical procedures when used in conjunction with a compatible RF/HF generator.

The joimax® bipolar electrosurgical instruments have not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures and should not be used for these procedures.

Model JVP32024 Vaporflex® Probe, Bipolar, Ball Tip 320

Model JVP27525S Vaporflex® Probe, Bipolar, Ball Tip 275

Model JVP25024 Vaporflex® Probe, Bipolar, Ball Tip 250

Model JBPP27020 Legato® Probe, Bipolar, Ball Tip 270

Model JBPH352505 Legato® Handpiece, Bipolar, for SurgiMax® / Endovapor® / Endovapor®²

Model JVH32030 Vaporflex® Handle and Shaft, 320

Model JVHS32030 Vaporflex® Shaft 320

Model JVHS27527S Vaporflex® Irrigation Shaft 275

Model JVHS25030 Vaporflex® Shaft 250

Model JVUH10010 Vaporflex® Handle

Model JVC35010 Vaporflex® Cable, for SurgiMax® / Endovapor®

Model JVC35020 Vaporflex® Cable, for Endovapor®²

Model JVHP10023 Vaporflex® Spare Sealing Rings

CONTRAINDICATIONS

- Do not use on patients for whom an electrosurgical procedure is contraindicated for any reason.
- Not for use in cardiac interventions.
- Do not use on patients with heart pacemakers or other electronic device implants without first consulting a qualified professional.


TECHNICAL DATA

- The maximum permitted output voltage is 1200 Vp at a crest factor of more than 1.9
- The maximum permitted output voltage is 650 Vp at a crest factor of less than 1.9

COMPONENTS

- (1) Vaporflex® Handle with (4) Sealing Ring, (5) Sealing Cap and (6) Locking Button
- (2) Shaft: available in different lengths and versions
Irrigation shaft: shaft version with Luer-Lock connector sideways to connect rinsing if needed
- (3) RF/HF-Cable 3.5 m long; the plug on cable is compatible with the following RF/HF generators:
joimax® (Endovapor®, Endovapor®2), **elliquence® (SurgiMax®, SurgiMax® plus, SurgiMax® Vapor)**

ASSEMBLY INSTRUCTIONS

- Ensure that the (4) Sealing Ring and the (5) Sealing Cap are securely attached to the (1) Vaporflex® Handle.
- Insert the bipolar probe (sterile disposable product) from a distal direction into the (1) handle until the markings of the probe are visible at the proximal end of the handle. The (6) Locking Button must be kept pressed.
- Release the (6) Locking Button and slowly draw the probe in a distal direction until it locks in place. Then check that it is securely fitted. Please mind the position of the probe as it can be locked in two directions (upwards or downwards).
- Guide the (2) shaft from a distal direction across the probe and firmly screw it to the (1) handle.
- Ensure that the RF/HF generator provides an ellman/elliquence® type bipolar 2-pin socket to connect the Vaporflex® bipolar instrument.
- Connect the (3) HF-cable to the (1) handle, the markings (blue dots) on the handle and cable must align.
- Connect the (3) HF-cable of the device to the RF/HF generator and only then switch on the RF/HF generator.

Before performing the procedure, the **Vaporflex® Handle** should be tested extracorporeal on a piece of moist gauze, for instance, to ensure that it is functioning. Upon activating the probe for a short duration, air bubbles should appear in the liquid on the tip of the probe.

DISASSEMBLY INSTRUCTIONS

- a. Switch off the RF/HF generator, remove the (3) RF/HF-cable first and then remove the plug from the handle.
- b. Unscrew the (2) shaft and after activating the (6) Locking Button, pull it in a distal direction from the handle together with the probe.
- c. Pull the probe (sterile disposable product) out from the (2) shaft and dispose it in accordance with internal or statutory regulations, since the probe is a biological hazard. Dispose the used RF/HF-probe as clinical waste in accordance with hospital guidelines.
- d. Refer to chapter cleaning, disinfection and sterilization for care of the **joimax® Vaporflex® Handle**.

Hold the RF/HF-cable only by the connector to plug in and out.

FUNCTION TEST

By squeezing the handle, the tip of the probe is pushed out of the shaft; the probe is activated only with the foot switch of the RF/HF generator.

After assembly or cleaning / disinfection, check all components for completeness, function, corrosion, damaged surfaces, chipping and cleanliness.

Take special care of a straight shaft without any sharp edges at the distal end. The seals on the handle are present and mounted. The cable is intact and free of kinks.

Inspect especially the cables for damage **prior to each use**, especially the insulations, visually under magnification or with a high voltage insulation testing device. Components with any defect are not allowed to be used anymore.

In case of a negative function test, reassemble the instrument. When proper function cannot be achieved, replace the faulty parts with new ones.

MATERIAL DURABILITY

When selecting cleaners and disinfectants, always ensure that they do not contain the following components:

- organic, mineral or oxidating acids (minimum allowed pH-value 5.5)
- strong lye (maximum allowed pH-value 11; we recommend a neutral / enzymatic or slightly alkaline cleaner)
- organic solvents (e.g. alcohol, ether, ketone, benzene)
- oxidants (e.g. hydrogen peroxides)
- halogens (chlorine, iodine, bromine)
- aromatic / halogenated hydrocarbons
- heavy metal salts

Never clean medical devices with metal brushes or steel wool!

PRODUCT LIFETIME

The lifetime of instruments is limited, because they are subject to a natural wear and aging.

Inappropriate handling e.g. overloading, unsuitable cleaning agents, high sterilization temperatures, etc. reduce the lifetime of the instruments.

A thorough examination of the instrument according to the function test must take place during the reprocessing process, after the disinfection and before the sterilization. If a deterioration / defect is determined, then the instrument is to be segregated immediately and may not further be used. The use of damaged and / or dirty instruments lies in the responsibility of the user. If you have questions contact joimax® GmbH or one of its representatives.

joimax® refuses any liability in the event of non-observance.

NOTES ON THE CORRECT USE OF THE DEVICE

- The Vaporflex® Handle will operate only with a bipolar probe of joimax® GmbH. Handle and probe must be assembled before use according to the instructions.
- Use the shaft only with a compatible probe. Note the different lengths.
- Select the cable according to the RF/HF generator available. Note the different plugs.
- Image interference on the monitor can be minimized if
 - the RF/HF cable is not led directly parallel with the camera cables,
 - the RF/HF cable is not looped,
 - the RF/HF generator is not switched on until the assembled Vaporflex® Handle has been connected.
- The power output of the RF/HF generator must be checked before each use and set up for the following application. Please follow the instructions for use of the RF/HF generator and the general instructions for electrosurgery.

WARNINGS AND POTENTIAL RISKS

- The use of non-flammable agents for cleaning and disinfection is recommended wherever possible; flammable agents for cleaning, disinfecting, or as solvents of adhesives shall be allowed to evaporate before application of RF/HF surgery; there is a risk of pooling of flammable solutions under the patient or in body depressions and in body cavities; fluids pooled in body depressions and cavities should be mopped up before RF/HF devices are used.
- Incorrect use of RF/HF current can cause endogenous and exogenous burns and explosions. Check before use that the power output set on the RF/HF generator is not higher than required for the surgical procedure. Details can be found in the instructions for use of the RF/HF generator.
- Before increasing power output because of insufficient coagulation, always check whether:
 - all connection contacts on the Vaporflex® Handle are in order (connector, handle, probe, shaft, etc.),
 - the foot switch is functioning perfectly,
 - the insulation on the RF/HF cable and the active probe is intact.
- joimax® probes may still be hot immediately after deactivation. Probes that have just been used, as well as temporarily unused probes, must not be placed on the patient or other flammable materials. They must be stored in a separate, insulated place.
- Only non-combustible gases (e.g. CO₂) can be used when using the bipolar probe during a procedure under insufflation.
- The functions that occur with proper use of the RF/HF generator can cause some materials, such as cotton and gauze, to ignite if they are saturated with oxygen.
- Connect the device to the RF/HF generator only when it is not activated.
- Do not activate the instrument when not in contact with target tissue.
- Ensure that insulated parts do not come in contact with hard, pointed or heavy objects during preparation, as such objects may damage the insulation and render the product unsuitable.
- The use of flammable anesthesia or gases such as laughing gas (N₂O) and oxygen should be avoided if an operation is being performed in the area of the thorax or head, unless these substances are sucked off.

CAUTIONS

- Read the instructions for use which are provided with the RF/HF-generator prior to using this device!
- Ensure that all package inserts including cautions and instructions for use are read and understood.
- Carefully read the instructions for use of the Vaporflex® Handle in full before surgical use.
- The components must be assembled in accordance with the specifications in the assembly instructions.

- The products mentioned above can be used only in medical facilities by a doctor familiar with their use and in accordance with the current medical standard.
- To avoid damaging the probe and the endoscope, the instrument may only be introduced into and withdrawn from the endoscope when the instrument is in a retracted position. In doing so, the handle should **not** be squeezed!
- To avoid damaging the probe and shaft, the instrument may only be activated in an extended position. In doing so, the handle must be squeezed!
- The power output of the RF/HF generator should always be set as low as necessary to achieve the desired effect.
- Avoid power output settings of the RF/HF generator where the maximum output voltage may exceed the maximum permitted voltage of the Vaporflex® Handle.
- Keep the active probes clean. Build-up of eschar may reduce the instrument's effectiveness. Do not activate the instrument while cleaning.

ELECTRO MAGNETIC COMPATIBILITY (EMC)

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents of the RF/HF generator. Portable and mobile RF/HF communication equipment can affect other medical electrical equipment.

The joimax® Vaporflex® Handle connected to a RF/HF generator [joimax® Endovapor®, joimax® Endovapor® 2, elliquence® (SurgiMax®, SurgiMax® plus, SurgiMax® Vapor)] is in compliance with the requirements of 6.1 and 6.2 of standard IEC 60601-1-2.

The use of accessory with medical electrical equipment (probe, handle, cable) other than those specified may result in increased emissions or decreased immunity of the medical electrical equipment.

PRINCIPLE OF BIPOLAR ELECTROSURGERY

Bipolar RF/HF surgery is the treatment of tissue by applying a high-frequency alternating current, which generates heat. Heat can shrink, cut, coagulate, and stop tissue bleeding. In bipolar configuration the RF/HF current flows from one branch of the instrument across the soft tissue to the second branch with a localized effect as required. The use of a neutral electrode is redundant with bipolar RF/HF applications. Branches are the non-insulated areas of the probe tip.

LIABILITY AND WARRANTY

This product is in accordance with MDD 93/42 EEC. Liability or warranty is excluded,

- if the notices and instructions in the instructions for use are not observed,
- if the product is not used according to regulations,
- if original parts are not used.

joimax® GmbH accepts no guarantee for products malfunctioning in connection with disinfectants and the cleaning and preparation procedure applied, including the effectiveness of the disinfectant.

LIMITED WARRANTY

joimax® guarantees to manufacture the products with utmost care.

joimax® GmbH issues the legally specified warranty period with regard to the function of instruments. The validity period of this warranty is limited to claims made within the specified warranty period after the instrument's purchase date, with reference to repairs, if required, and with specification of the invoice number.

This warranty only encompasses defects not caused by normal wear and tear, misuse, wrong handling, no or incorrect preparation or force majeure. This warranty excludes wearing parts. This is the only valid warranty and replaces all other issued warranty declarations.

joimax® component sets as well as individual products are compatible with each other. Before using individual joimax® products / sets with other products, the user must ensure the user-specific compatibility of individual products. We recommend introductory training and regular attendance at training events. We are available for any further required information. joimax® has no influence on the application of the product, on the diagnosis of the patient and handling of the product outside of joimax®. joimax® cannot guarantee good results nor the complication free application of the products. joimax® employees are not entitled to change the previously mentioned conditions, extend liability or enter additional product related obligations. joimax® reserves the right to product changes. This document contains information protected by copyright and property law and may not be copied in full or in parts thereof or transferred to a further medium in any form, without express written permission from joimax® GmbH. Distribution to third parties is prohibited. The products and names used here may be the registered brands of joimax® or other companies.

PRIOR TO USE

- The Vaporflex® Handle is delivered non-sterile and must be cleaned and sterilized before its first surgical use. Carefully read these instructions before surgical use and please note the specifications of joimax® GmbH.
- A complete sterile joimax® Vaporflex® Handle with probe should always be available as substitute.

A visual and functional inspection must be carried out prior to every use!

INSTRUCTIONS FOR USE

1. Assemble joimax® Vaporflex® Handle according to the instructions provided.
2. Select the required program on the HF/RF generator.
3. Adjust the power setting.

It is recommended to start the RF/HF generator with low energy settings. If required, the energy can be gradually increased. Progress of work must be regularly monitored.

4. Perform a function test as described.
5. Insert the joimax® Vaporflex® instrument into the appropriate endoscopic working channel.
6. Once the probe tip contacts the tissue to be treated, press the foot switch to activate the connected RF/HF generator.

PROCESSING INSTRUCTIONS

PRETREATMENT

Soaking must be carried out immediately after use (within a maximum of 2 hours). Use water or an aldehyde-free disinfectant for this purpose (to avoid fixation of blood contamination) with proven effectiveness (FDA cleared), which are suitable for disinfecting instruments and which are compatible with the instruments.

A synthetic fleece or plastic cleaning brush can be used if necessary for the preliminary removal of residue.

Do not pre-clean RF/HF cables and hand pieces in the ultrasonic bath.

A disinfectant used for soaking is only intended for protecting personnel and cannot replace the later disinfecting step after cleaning has been carried out.

MANUAL CLEANING AND DISINFECTION – FDA CLEARED

CLEANING

The following instructions only apply to the Vaporflex® Handle, and not for the Vaporflex® Bipolar Probe, which is single-use.

1. Disassemble the device completely following the disassembly instructions.
2. Submerge the disassembled parts into the cleaning bath for the specified reaction time, so that the instruments are sufficiently covered (careful brushing with a soft brush, if required), ensuring that the instruments do not touch each other.

The time that devices are soaked in the cleaning bath is given by the manufacturer of the cleaning agent (for example, the manufacturer of Cidezyme® recommends a minimum of 1 minute).

3. At the beginning and at the end of the reaction time, rinse all luminae at least five times by using a disposable syringe (minimum volume 100 ml) and the Luer-Lock connection.
4. Then remove the medical devices from the cleaning bath and thoroughly rinse them with water at least three times.
5. Ensure that the device is thoroughly cleaned and visually clean.
6. Repeat (2.) if necessary.

DISINFECTION

7. Submerge the disassembled, cleaned and checked medical devices into the disinfection bath for the specified reaction time, so that the medical devices are sufficiently covered, ensuring that the instruments do not touch one another. Insert the cable loose, it may not be kinked or crushed.

Limit the time that devices are soaked in the disinfectant solution to the minimum time recommended by the disinfectant manufacturer (for example, the manufacturer of Cidex® OPA recommends a minimum of 12 minutes).

8. At the beginning and at the end of the reaction time, rinse all instrument luminae five times by using a disposable syringe (minimum volume 100 ml) and the Luer-Lock connection or by applying a rinsing hose.
9. Remove the instruments from the disinfection bath and thoroughly rinse them with purified water at least five times.
10. Dry the instruments by blowing off and out with filtered pressurized air.

Only dry the products with compressed air at pressures below 3 bar to avoid possible damage.

11. Pack the instruments, preferably immediately after removal (refer to packing), after additional drying in a clean place.

Components which are still soiled after the reprocessing procedure, must be cleaned and disinfected again.

joimax® recommendation regarding cleaning and disinfection:

Proof of fundamental instrument suitability for effective manual cleaning and disinfection was provided by an independent accredited testing laboratory with the use of the cleaner Cidezyme® and the disinfectant CIDEX® OPA, whereby the procedures described above were used.

AUTOMATED CLEANING AND DISINFECTION – NOT FDA CLEARED, FOR OUTSIDE US MARKETS ONLY

When selecting the disinfectant, ensure the following

- that the disinfectant has validated efficacy (e.g. DGHM or FDA approval or CE Mark according to DIN EN ISO 15883),
- that, if possible, a validated thermal disinfection program (A0 value > 3000 or, or, in older devices, at least 5 min at 90°C) is used
- that the used program is suitable for the instruments and has sufficient rinsing cycles,
- that only sterile or germ-reduced, (max. 10 germs / ml) as well as, endotoxine reduced (max. 0.25 endotoxine units / ml) water (e.g. purified / highly purified water) is used for rinsing.
- that the air used for drying is filtered
- that the disinfectant is maintained and inspected regularly.

Caution: Chemical disinfection poses the hazard of disinfectant residues on the instruments.

When selecting the applied cleaning system, ensure

- that it is suitable for cleaning metal and plastic instruments,
- that – unless thermal disinfection is used – a suitable disinfectant with validated efficacy (e.g. VAH / DGHM or FDA approved or with a CE Mark) is used and that this is compatible with the applied cleanser and
- that the applied chemicals are compatible with the instruments (refer to chapter 12 on „Material durability“).

Caution: The concentration provided by the manufacturer of the cleaning resp. disinfection agent must be followed.

Procedure:

1. Disassemble the medical devices completely.
2. Place the disassembled medical devices into the disinfectant, ensure that the instruments do not touch each other. Insert the cable loose, they may not be kinked or crushed.
If applicable: Connect all instrument lumens to the rinsing connection by using the Luer-Lock connection on the product or by applying a suitable rinsing hose to the rinsing adapter of the disinfectant. Otherwise, rinse all lumens of the instruments five times using a disposable syringe (minimum volume 10-50 ml, depending on the instrument's size) slipped over the Luer connector or by a tube- at inlet and outlet of the lumen.
3. Start the program.
4. When the program is finished, remove the medical devices from the disinfectant
5. Inspect and pack the medical devices preferably immediately after removal after additional drying in a clean place.

joimax® recommendation regarding cleaning and disinfection: Proof of fundamental instrument suitability for effective mechanical cleaning and disinfection was provided by an independent accredited testing laboratory with the use of the disinfectant G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent Neodisher mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg), whereby the procedure described above was used.

PACKING

Sort the cleaned and disinfected instruments into the provided sterilization tray. Then pack the sterilization tray into disposable dual packaging and / or sterilization containers complying with the following requirements:

- DIN EN ISO / ANSI AAMI ISO 11607-1; 11607-2
- suitable for steam sterilization (heat resistant minimum 141° C / 286° F, sufficient vapor permeability)
- regularly maintained
- cleared by the FDA (No. K803293, K953776 and K973827)

The packing is to be marked with information about method and date of sterilization, expiration date and batch.

STERILIZATION

The individual user is responsible for the validation of the specified sterilization process. The following method must be used: Steam sterilization with fractionated vacuum procedure with a temperature of 132°C/270°F for a minimum of 3 minutes / a maximum of 5 minutes and a drying time of 20 minutes.

Specific sterilization parameters for a particular unit should be validated by the use of spore strips or inoculated products. Any manufacturer's sterilization instructions for a particular unit should take precedence over recommendations given here.

- Do not exceed 138°C/280°F
- Do not use hot air sterilization
- The FDA does not recognize sterilization methods that have not been 510(k)-cleared

joimax® recommendation regarding sterilization:

Proof of basic instrument suitability for effective steam sterilization was performed by an independent accredited testing laboratory by using the steam sterilizer Getinge 400HC and 500HC (FDA clearance K103504) with the fractional vacuum procedure. Typical conditions prevalent at a clinic or physician's practice as well as the above-mentioned procedure were taken into consideration.

STORAGE

Temperature: -20° C to +50° C / -4° F to +122° F

Atmospheric pressure: 50 –106 kPa

Relative humidity: 0 – 75 %, non-condensing

CUSTOMER SERVICE

In case of questions or problems please contact the customer service of joimax® GmbH.

In the interest of rapid processing of enquiries please send us the product with the following information:

- Article number (REF)
- Serial number (SN or LOT)
- Description of the error as accurate as possible


Please notify us of any events according to MPG § 29 or § 2 of the Medical Devices Safety Plan Ordinance under the following telephone number:

for the US 949-859-3472 (24-hour availability)

for all other countries +49 (0) 721 255 14-999 (24-hour availability).

The CE Mark confirms product compliance with the Medical Device Directive 93/42/EEC and the German Medical Devices Act (MPG).

SIGNS AND SYMBOLS

	Caution
	Follow instructions for use
	Manufacturer
	Batch code
	Catalogue number
	CE marking in accordance with Directive MDD 93/42/EEC
	Caution: U.S. Federal Law restricts this device to sale by, or on the order of, a physician
	Quantity one each
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation

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joimax®, Inc.