

Versicon[®]
Versatile Irrigation Control



**FOR ENDOSCOPIC SPINAL SURGERY
AND ARTHROSCOPY**



**THE UNIQUE
SPINE PUMP!**



**IRRIGATION PUMP
WITH MULTIPLE FUNCTIONS**

- Integrated spine mode for low flow and pressure
- Permanent control of flow and pressure
- Rapidly insertable, disposable tubing set
- Replaceable patient line with check valve





TECHNICAL SPECIFICATIONS

SPECIFICATIONS SPINE PORGRAM

Flow Range	0 – 750 ml/min
Pressure Range	10 – 100 mmHG

SPECIFICATIONS ARTHRO PORGRAM

Flow Range	0 – 2500 ml/min
Pressure Range	10 – 200 mmHG


GENERAL SPECIFICATIONS

Temperature Range	10°C – 35°C (50°F - 95°F)
Humidity Range	30 – 70 %
Power Supply	100 – 230 V / AC, 50 – 60 Hz
Main Fuse	2 x T3.15 AL 250 V (5 x 20 mm)
Power Consumption	60 VA
Dimensions in mm	(H x W D) 95 x 380 x 320
Weight	10 kg (22lb)
Safety Class	IEC 60601-1
Water Protection Class	IPX 1

ORDER INFORMATION

REF	DESCRIPTION
JISP3000	joimax® Versicon® spine pump incl. power cable, replacement fuses, sensor protection and instruction for use
JTSB350D	joimax® Versicon® tubing set (single-use)
JPPL0350	joimax® Versicon® patient line, 3.5 m (single-use)
FSSC00001	joimax® stopcock for Foraminoscope®
SCLL1520	joimax® Versicon® sealing cap
JIFS0001	joimax® Versicon® footswitch (optional)

DISPOSABLE TUBING SET



joined minimal access

joimax® GmbH
Amalienbadstrasse 41, RaumFabrik 61
76227 Karlsruhe, Germany

Phone +49 (0) 721 255 14-0
Fax +49 (0) 721 255 14-920
E-Mail info@joimax.com
Net www.joimax.com

joimax®, Inc.
140 Technology Drive, Suite 150
Irvine, CA 92618, USA

Phone +1 949 859 3472
Fax +1 949 859 3473
E-Mail info@joimaxusa.com
Net www.joimax.com

joimax® Asia
Rykadan Capital Tower,
135 Hoi Bun Road,
Kwun Tong, Hong Kong

Phone +852 29116418
E-Mail asia@joimax.com
Net www.joimax.com

joimax® ASEAN
Regus Vision Exchange
2 Venture Drive, Level 24-01 - 24-32
Singapore 608526

Phone +65 6914 9227
E-Mail asia@joimax.com
Net www.joimax.com